Split Stream® Catheter
 Καθετήρ Split Stream®
 Kateter Split Stream®
 Katetr Split Stream®
 Split Stream®-kateter
 Split Stream® katheter
 Split Stream® kateeter
 Split Stream® -katetri
 Cathéter Split Stream®
 Split Stream® Katheter
 Καθετήρας Split Stream®
 Split Stream® katéter
 Catetere Split Stream®
 Split Stream® katetrs
 „Split Stream®“ kateteris
 Cewnik Split Stream®
 Cateter Split Stream®
 Cateter Split Stream®
 Cathéter Split Stream®
 Kateter Split Stream®
 Cathéter Split Stream®
 Split Stream®-kateter
 Split Stream® Kateter

TABLE OF CONTENTS

en English	1	it Italian	25
bg Bulgarian	3	lv Latvian	27
hr Croatian	5	lt Lithuanian	29
cz Czech	7	pl Polish	31
da Danish	9	pt Portuguese	33
nl Dutch	11	ro Romanian	35
et Estonian	13	sk Slovak	37
fi Finnish	15	sl Slovenian	39
fr French	17	es Spanish	41
de German	19	sv Swedish	43
el Greek	21	tr Turkish	45
hu Hungarian	23		

en English

Split Stream® Catheter Patient Information Packet

PIP Statement:

This patient information packet is intended to provide general information regarding this device and should not replace the Information for Use.

Information for Users/Healthcare Professionals:

The following information is intended for users/healthcare professionals. Following this information, there is a summary intended for patients.

Indications for Use:

The Split Stream® Catheter is indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis.

Device Lifetime:

Medcomp® catheters are subjected to, and must pass, simulated use testing intended to replicate use 3 times per week for 12 months as part of device development. The Split Stream® Catheter passed this testing. Although Medcomp® catheters contain no materials which degrade over time, fully functional catheters may be removed

for other reasons, such as intractable infection, change of therapy (such as Renal replacement (transplant) or use of an arterio-venous graft/fistula). Published clinical literature does not always focus on the physical lifetime of a catheter for these reasons. In the case of the Split Stream® Catheter, 45, 49, and 47 catheters had 112, 130.7, and 130.5 day mean durations of use that has been found in clinical use reported to date¹.

Based on this information, the Split Stream® Catheter has a 12-month lifetime; however, the decision to remove and/or replace the catheter should be based on clinical performance and need, and not any predetermined point in time.

MRI Safety Information:

MR Conditional – 1.5 Tesla and 3 Tesla (artifacts may present imaging problems if MRI area of interest is on or near area where device is located)

Report Conclusion: MRI Information MR Conditional

Non-clinical testing has demonstrated that the Split Stream Catheter is MR Conditional. A patient with this device can be safely scanned in an MR system with the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0-Tesla (3.0T).
- Maximum spatial gradient field of 19 T/m (1900 G/cm).
- Maximum MR system reported, whole-body averaged specific absorption rate (SAR) of 2.0 J/kg (normal operating mode).

Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than or equal to 2.5°C after 15 minutes of continuous scanning. In non-clinical testing, the image artifact caused by the device extends up to 5.5 cm when imaged with a gradient echo pulse sequence in a 3.0T MR system.

Site Care:

Clean skin around catheter. Chlorhexidine gluconate solutions are recommended; however, iodine-based solutions can also be used. Cover the exit site with occlusive dressing and leave extensions, clamps, and caps exposed for access by staff.

Wound dressings must be kept clean and dry. **Patients must not swim, shower, or soak dressing while bathing.** If profuse perspiration or accidental wetting compromises adhesion of dressing, the medical or nursing staff must change the dressing under sterile conditions.

Alternate compatible solutions/ointments include:

- ChlorPrep®, Chlorhexidine Gluconate 2%
- Betadine®, Povidone Iodine 10%
- Isopropanol 70%
- Hydrogen Peroxide 3%
- Anasept®, 0.057% Sodium Hypochlorite
- Bactroban®, Mupirocin Ointment 2%
- Silvadene®, 1% Silver Sulfadiazine Cream

Solutions should be allowed to completely dry before applying an occlusive dressing.

Warnings, Precautions, or Measures to Be Taken by Healthcare Professional:

- Care must be taken when using sharp objects or needles in close proximity to the catheter lumen. Contact from sharp objects may cause catheter failure.
- Do not clamp the dual lumen portion of the catheter. Clamp only the extensions. Do not use serrated forceps; use only the in-line clamps provided.
- Extension clamps should only be open for aspiration, flushing, and dialysis treatment.
- Always review hospital or unit protocol, potential complications and their treatment, warnings, and precautions prior to undertaking any type of mechanical or chemical intervention in response to catheter performance problems.
- Due to the risk of exposure to HIV (Human immunodeficiency Virus) or other blood borne pathogens, health care professionals should always use Universal Blood and Body Fluid precautions in the care of all patients. To prevent accidents, assure the security of all caps and bloodline connections prior to and between treatments.
- Call 215-256-4201 or visit <http://www.medcomp.net> for any information to ensure safe use of the device, warnings, precautions or measures to be taken by the patient or a healthcare professional.

How To Fill Out Patient ID Card:

Instructions for Completion of Patient ID Card:

1. Name of the patient or patient ID. To be filled by the healthcare institution/provider.
2. Date of implantation. To be filled by the healthcare institution/provider.
3. Name and address of the healthcare institution/provider. To be filled by the healthcare institution/provider.

www.medcomp.net/patientinformation

Name: _____

Date of Implantation: _____

Name and Address of Healthcare Institution/Provider: _____

Medical Components, Inc.
 1997 Park Dr., Hollidown, PA 19108, U.S.A.
 P: 215-256-4201 F: 215-256-1787
www.medcomp.net

PATIENT INFORMATION:

The information presented below is intended for patients or lay persons. A more extensive summary of this device prepared for healthcare professionals is found in the first part of this document.

Contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

Site Care:

Clean skin around catheter. Chlorhexidine gluconate solutions are recommended; however, iodine-based solutions can also be used. Cover the exit site with occlusive dressing and leave extensions, clamps, and caps exposed for access by staff. Wound dressings must be kept clean and dry. **Patients must not swim, shower, or soak dressing while bathing.** If profuse perspiration or accidental wetting compromises adhesion of dressing, the medical or nursing staff must change the dressing under sterile conditions.

Alternate compatible solutions/ointments include:

- ChlorPrep®, Chlorhexidine Gluconate 2%
- Betadine®, Povidone Iodine 10%
- Isopropanol 70%
- Hydrogen Peroxide 3%
- Anasept®, 0.057% Sodium Hypochlorite
- Bactroban®, Mupirocin Ointment 2%
- Silvadene®, 1% Silver Sulfadiazine Cream

Solutions should be allowed to completely dry before applying an occlusive dressing.

Warnings, Precautions, or Measures to Be Taken by Patient:

- To reduce the risk of bacteria entering the catheter, wear a mask over your nose and mouth whenever the catheter is accessed.
- Keep the catheter dressing clean and dry. The dressing should be changed by a medical professional at each dialysis session.
- Avoid letting the catheter or catheter site go under water. Moisture near the catheter site can potentially lead to an infection.
- Ask the doctor to explain the signs and symptoms of catheter infection.
- Never remove the cap at the end of the catheter. The cap and clamps of the catheter shall be kept closed when not being used for dialysis.

www.medcomp.net/patientinformation

Materials or Substances the Patient May Be Exposed to:

The percentage ranges in the table below are based on the weight of the assembled 14F x 24 cm (20.32 g) and 14F x 40 cm (22.98 g) Split Stream catheters.

14F Split Stream®	
Material	% Weight (w/w)
Polyurethane	56.62 – 59.41
Acetal copolymer	15.01 – 17.33
Polypropylene	10.45 – 11.07
Silicone	5.31 – 6.13
Barium sulfate	6.31 – 8.14
Stainless Steel	2.65 – 3.06
Nylon	4.30 – 4.96
Acrylonitrile Butadiene Styrene	3.52 – 4.07
Polyethylene terephthalate	1.65 – 1.91

The percentage ranges in the table below are based on the weight of the assembled 16F x 24 cm (20.13 g) and 16F x 40 cm (24.16 g) Split Stream catheters.








16F Split Stream®	
Material	% Weight (w/w)
Polyurethane	49.86 – 53.37
Acetal copolymer	10.37 – 11.74
Polypropylene	9.24 – 10.45
Silicone	5.01 – 5.67
Barium sulfate	5.69 – 7.36
Stainless Steel	2.37 – 2.68
Nylon	8.67 – 7.66
Acrylonitrile Butadiene Styrene	3.14 – 3.56
Polyethylene terephthalate	1.48 – 1.68

Note: Accessories containing stainless steel may contain up to 0.4% weight of the CMR substance cobalt.

References:

- Silva J, e Costa T, Baptista A, Ramos A, Ponce P. Catheter-Related Bacteremia in Hemodialysis: Which Preventive Measures to Take? Nephron Clin Pract. 2008;110(4):c251-c257.

SYMBOL TABLE

	Manufacturer*
	Serial Number*
	Batch Code*
	Authorized Representative in the European Community/European Union*
	Medical Device*
	Unique Device Identifier*
	MR Conditional - 1.5, 3.0 Tesla***
	Patient Information Website*
	Date*
	Healthcare Center or Doctor*
	Patient Identification*
	Date of Manufacture*
	United Kingdom Responsible Person**
	Contains Hazardous Substances*
	Contains the CMR Substance, Cobalt

*This symbol is in accordance with ISO 15223-1.

**Not a recognized symbol.

***This Symbol is in accordance with ASTM F 2503-20.

Medcomp® and Split Stream® are trademarks of Medical Components, Inc. registered in the United States.

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Harleysville, PA 19438, U.S.A.
MPS Medical-Product-Service UK Ltd
MPS Medical Product Service GmbH