6.5F Tesio® Catheter
Patient Information Packet

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6.5 Tesio® Catheter Patient Information Packet

PIP Statement:

This patient information packet is intended to provide general information regarding this device and should not replace the Information for Use.

Information for Users/Healthcare Professionals:

The following information is intended for users/healthcare professionals. Following this information, there is a summary intended for patients.

Indications for Use:

- The 6.5F Tesio® Catheter is indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis.

Device Lifetime:

Medcomp® catheters are subjected to, and must pass, simulated use testing intended to replicate use 3 times per week for 12 months as part of device development. The 6.5F Tesio® Catheter passed this testing. Although Medcomp® catheters contain no materials which degrade over time, fully functional catheters may be removed for other reasons, such as intractable infection, change

of therapy (such as Renal replacement (transplant) or use of an arterio-venous graft/fistula). Published clinical literature does not always focus on the physical lifetime of a catheter for these reasons. In the case of the 6.5F Tesio® Catheter, 11 catheters had a median 244 day [Range: 1 - 344 days] duration of use that has been found in clinical use reported to date*.

Based on this information, the 6.5F Tesio® Catheter has a 12-month lifetime; however, the decision to remove and/or replace the catheter should be based on clinical performance and need, and not any predetermined point in time.

MR! Safety Information:

 MR Conditional – 1.5 Tesla and 3 Tesla (artifacts may present imaging problems if MRI area of interest is on or near area where device is located)

Report Conclusion: MRI Information MR Conditional

The Tesio® Catheter was determined to be MR-conditional. Non-clinical testing demonstrated that the Implantable Vascular Access Port is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MR!-Related Heating

In non-clinical testing, the Tesio® Catheter (polyurethane with embedded stainless steel connector) produced the following temperature rise during MRI performed for 15-min in the 3- Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system:

Highest temperature change +1.6°C

Therefore, the MR!-related heating experiments for the Tesio® Catheter (polyurethane with embedded stainless steel connector) at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 -W/kg (i.e., associated with a calorimetry measured whole body averaged value of

2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.6°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Tesio® Catheter (polyurethane with embedded stainless steel connector). Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	778-mm²	233-mm²	1,456-mm²	1,778-mm²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Site Care:

Clean skin around catheter. Chlorhexidine gluconate solutions are recommended; however, iodine-based solutions can also be used. Cover the exit site with occlusive dressing and leave extensions, clamps, and caps exposed for access by staff. Wound dressings must be kept clean and dry. **Patients must not swim, shower, or soak dressing while bathing.** If profuse perspiration or accidental wetting compromises adhesion of dressing, the medical or nursing staff must change the dressing under sterile conditions.

Alternate compatible solutions/ointments include:

- ChloraPrep®, Chlorhexidine Gluconate 2%
- Betadine®, Povidone Iodine 10%
- Isopropanol 70%
- Hydrogen Peroxide 3%
- Anasept®, 0.057% Sodium Hypochlorite
- Bactroban®, Mupirocin Ointment 2%
- Silvadene®, 1% Silver Sulfadiazine Cream
- Providine, 10% Povidone - Iodine Ointment
- Polysporin®, Bacitracin Zinc, Polymyxin B Sulfate Ointment
- Gentamycin Sulfate Cream 0.1%

Solutions should be allowed to completely dry before applying an occlusive dressing.

Warnings, Precautions, or Measures to Be Taken by Healthcare Professional:

- Care must be taken when using sharp objects or needles in close proximity to


the catheter lumen. Contact from sharp objects may cause catheter failure.

- Extension clamps should only be open for aspiration, flushing, and dialysis treatment.
- Always review hospital or unit protocol, potential complications and their treatment, warnings, and precautions prior to undertaking any type of mechanical or chemical intervention in response to catheter performance problems.
- Due to the risk of exposure to HIV (Human immunodeficiency Virus) or other blood borne pathogens, health care professionals should always use Universal Blood and Body Fluid precautions in the care of all patients. To prevent accidents, assure the security of all caps and bloodline connections prior to and between treatments.
- Call 215-256-4201 or visit <http://www.medcomp.net> for any information to ensure safe use of the device, warnings, precautions or measures to be taken by the patient or a healthcare professional.

How To Fill Out Patient ID Card:

Instructions for Completion of Patient ID Card:

1. Name of the patient or patient ID. To be filled by the healthcare institution/provider.
2. Date of implantation. To be filled by the healthcare institution/provider.
3. Name and address of the healthcare institution/provider. To be filled by the healthcare institution/provider.

PATIENT INFORMATION:

The information presented below is intended for patients or lay persons. A more extensive summary of this device prepared for healthcare professionals is found in the first part of this document.

Contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

Site Care:

Clean skin around catheter. Chlorhexidine gluconate solutions are recommended; however, iodine-based solutions can also be used. Cover the exit site with occlusive dressing and leave extensions, clamps, and caps exposed for access by staff. Wound dressings must be kept clean and dry. **Patients must not swim, shower, or soak dressing while bathing.** If profuse perspiration or accidental wetting compromises adhesion of dressing, the medical or nursing staff must change the dressing under sterile conditions.

Alternate compatible solutions/ointments include:

- ChloraPrep®, Chlorhexidine Gluconate 2%
- Betadine®, Povidone Iodine 10%
- Isopropanol 70%
- Hydrogen Peroxide 3%
- Anasept®, 0.057% Sodium Hypochlorite
- Bactroban®, Mupirocin Ointment 2%
- Silvadene®, 1% Silver Sulfadiazine Cream
- Providine, 10% Povidone - Iodine Ointment
- Polysporin®, Bacitracin Zinc, Polymyxin B Sulfate Ointment
- Gentamycin Sulfate Cream 0.1%

Solutions should be allowed to completely dry before applying an occlusive dressing.

Warnings, Precautions, or Measures to Be Taken by Patient:

- To reduce the risk of bacteria entering the catheter, wear a mask over your nose and mouth whenever the catheter is accessed.
- Keep the catheter dressing clean and dry. The dressing should be changed by a medical professional at each dialysis session.
- Avoid letting the catheter or catheter site go under water. Moisture near the catheter site can potentially lead to an infection.
- Ask the doctor to explain the signs and symptoms of catheter infection.
- Never remove the cap at the end of the catheter. The cap and clamps of the catheter must be kept closed when not being used for dialysis.

www.medcomp.net/patientinformation

Materials or Substances the Patient May be Exposed to:

The percentages in the table below are based on the weight of the assembled catheters, 9.47g.







Material	% Weight (w/w)
Polyvinyl chloride	32.45
Acetal copolymer	40.63
Polyurethane	20.84
Barium sulfate	3.18
Polyethylene terephthalate	1.60
Silicone	1.30

Note: Accessories containing stainless steel may contain up to 0.4% weight of the CMR substance cobalt.

Reference:

- Sheth RD, Kale AS, Brewer ED, Brandt ML, Nuchtern JG, Goldstein SL. Successful use of Tesio catheters in pediatric patients receiving chronic hemodialysis. Am J Kidney Dis. 2001; 38(3):553-9.

SYMBOL TABLE

	Manufacturer *
	Serial Number *
	Batch Code *
	Authorized Representative in the European Community/European Union *
	Medical Device *
	Unique Device Identifier *
	MR Conditional - 1.5, 3.0 Tesla ****
	Patient Information Website *
	Date *
	Healthcare Center or Doctor *
	Patient Identification *
	Date of Manufacture *
	United Kingdom Responsible Person **
	Contains Hazardous Substances *
	Contains the CMR Substance, Cobalt

**This symbol is in accordance with ISO 15223-1.

***Not a recognized symbol.

****This Symbol is in accordance with ASTM F 2903-20.

Medcomp® and Tesio® are trademarks of Medical Components, Inc. registered in the United States.

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Medcomp ve Tesio® Medical Components, Inc. firmasının Amerika Birleşik Devletleri'nde kayıtlı ticari markalarıdır.
Harleysville, PA 19438, U.S.A.
MPS Medical-Product-Service UK Ltd
MPS Medical Product Service GmbH