Symetrex® Catheter
Patient Information Packet

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Symetrex® Catheter Patient Information Packet

PIP Statement:

This patient information packet is intended to provide general information regarding this device and should not replace the Information for Use.

Information for Users/Healthcare Professionals:

The following information is intended for users/healthcare professionals. Following that information is a summary intended for patients.

Indications for Use:

The Symetrex® Catheter is indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis.

Device Lifetime:

Medcomp® catheters are subjected to, and must pass, simulated use testing intended to replicate 12 months use as part of device development. The Symetrex® Catheter passed this testing. Although Medcomp® catheters contain no materials which degrade over time, fully functional catheters may be removed for other reasons, such as intractable infection or change of therapy (such as Renal replacement (transplant)

or use of an arterio-venous graft/fistula). Published clinical literature does not always focus on the physical lifetime of a catheter for these reasons. In the case of the Symetrex® Catheter, 46 catheters had an 87 day [95%CI: 61.7 – 112.3 days] duration of use that has been found in clinical use reported to date⁴.

Based on this information, the Symetrex® Catheter has a 12-month lifetime; however, the decision to remove and/or replace the catheter should be based on clinical performance and need, and not any predetermined point in time.

MRI Safety Information:

The Symetrex® Catheter is
MR-Conditional.

A patient utilizing this device can be safely scanned in an MR system with the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0-Tesla (3.0T)
- Maximum spatial gradient field of 19 T/m (1900 G/cm)
- Maximum MR system reported, whole-body averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)

Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than or equal to 2°C after 15 minutes of continuous scanning in non-clinical testing. The image artifact caused by the device extends radially up to 0.2cm from the device when imaged with a gradient echo pulse sequence in a 1.5T MR system or a spin echo pulse in a 3.0T MR system.

Site Care:

Clean skin around catheter. Chlorhexidine gluconate solutions are recommended; however, iodine-based solutions can also be used. Cover the exit site with occlusive dressing and leave extensions, clamps, and caps exposed for access by staff. Wound dressings must be kept clean and dry.

Patients must not swim, shower, or soak dressing while bathing. If profuse perspiration or accidental wetting compromises adhesion of dressing, the medical or nursing staff must change the dressing under sterile conditions.

Alternate compatible solutions/oointments include:

- ChloraPrep®, Chlorhexidine Gluconate 2%
- Betadine®, Povidone Iodine 10%
- Isopropanol 70%
- Hydrogen Peroxide 3%
- Anasept®, 0.057% Sodium Hypochlorite
- Bactroban®, Mupirocin Ointment 2%
- Silvadene®, 1% Silver Sulfadiazine Cream
- Providine, 10% Povidone – Iodine Ointment
- Polysporin®, Bacitracin Zinc, Polymyxin B Sulfate Ointment
- Gentamycin Sulfate Cream 0.1%

Solutions should be allowed to completely dry before applying an occlusive dressing.

Warnings, Precautions, or Measures to be Taken by Healthcare Professional:

- Care must be taken when using sharp objects or needles in close proximity to the catheter lumen. Contact from sharp objects may cause catheter failure.
- Do not clamp the dual lumen portion of the catheter. Clamp only the extensions. Do not use serrated forceps; use only the in-line clamps provided.
- Extension clamps should only be open for aspiration, flushing, and dialysis treatment.
- Always review hospital or unit protocol, potential complications and their treatment, warnings, and precautions prior to undertaking any type of mechanical or chemical intervention in response to catheter performance problems.
- Due to the risk of exposure to HIV (Human immunodeficiency Virus) or other blood borne pathogens, health care professionals should always use Universal Blood and Body Fluid precautions in the care of all patients.
- Call 215-256-4201 or visit www.medcomp.net for any information to ensure safe use of the device, warnings, precautions or measures to be taken by the patient or a healthcare professional.
- The CMR substance Cobalt is a naturally occurring component of stainless steel. Based on biocompatibility evaluation it was determined that the main hazards of stainless steels are related to the processing of the material, especially welding, thus not applicable to the intended use of the device. Stainless steels used in these devices are unlikely to reach exposure levels that will elicit carcinogenicity, mutagenicity or reproductive toxicity.

How To Fill Out Patient ID Card:

Instructions for Completion of Patient ID Card:

1. Name of the patient or patient ID. To be filled by the healthcare institution/provider.
2. Date of implantation. To be filled by the healthcare institution/provider.
3. Name and address of the healthcare institution/provider. To be filled by the healthcare institution/provider.

The form is titled "Patient ID Card" and includes a URL <http://www.medcomp.net/patientinformation>. It has three main sections for text entry: "Name", "Date", and "Address". Below these is a logo for "Medical Components, Inc." with contact information: "(800) 646-1234, (215) 256-4201" and "www.medcomp.net".

PATIENT INFORMATION:

The information presented below is intended for patients or lay persons. A more extensive summary of this device prepared for healthcare professionals is found in the first part of this document.

Contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

Site Care:

Clean skin around catheter. Chlorhexidine gluconate solutions are recommended; however, iodine-based solutions can also be used. Cover the exit site with occlusive dressing and leave extensions, clamps, and caps exposed for access by staff. Wound dressings must be kept clean and dry.

Patients must not swim, shower, or soak dressing while bathing. If profuse perspiration or accidental wetting compromises adhesion of dressing, the medical or nursing staff must change the dressing under sterile conditions.

Alternate compatible solutions/ointments include:

- ChloraPrep®, Chlorhexidine Gluconate 2%
- Betadine®, Povidone Iodine 10%
- Isopropanol 70%
- Hydrogen Peroxide 3%
- Anasept®, 0.057% Sodium Hypochlorite

- Bactroban®, Mupirocin Ointment 2%
- Silvadene®, 1% Silver Sulfadiazine Cream
- Providine, 10% Povidone – Iodine Ointment
- Polysporin®, Bacitracin Zinc, Polymyxin B Sulfate Ointment
- Gentamycin Sulfate Cream 0.1%

Solutions should be allowed to completely dry before applying an occlusive dressing.

Warnings, Precautions, or Measures to be Taken by Patient:

- It's important to understand what blood flow rate, Kt/V, and urea reduction ratio (URR) mean to you and your dialysis session. Ask a clinician what these mean and what numbers are normal for you.
- If a decrease in blood flow rate occurs during more than one treatment in a week, contact your doctor; the catheter should be checked.
- To reduce the risk of bacteria entering the catheter, wear a mask over your nose and mouth whenever the catheter is accessed.
- Keep the catheter dressing clean and dry. The dressing should be changed by a medical professional at each dialysis session.
- Avoid letting the catheter or catheter site go under water. Moisture near the catheter site can potentially lead to an infection.
- Ask the doctor to explain the signs and symptoms of catheter infection. Never remove the cap at the end of the catheter. The cap and clamps of the catheter shall be kept closed when not being used for dialysis.

www.medcomp.net/patientinformation

Materials or Substances the Patient May be Exposed to:

The percentage ranges in the table below are based on the weights of the 15.5F x 19cm Double Lumen catheter without sideholes (11.88g), 15.5F x 19cm Double Lumen catheter with sideholes (11.87g), 15.5F x 42cm Double Lumen catheter without sideholes (15.60g), and 15.5F x 42cm Double Lumen catheter with sideholes (15.59g).

Material	% Weight (w/w)
Polyurethane	35.96 - 46.50
Acetyl Copolymer	15.82 - 20.79
Isoplast	20.94 - 27.53
Barium Sulfate	7.69 - 10.64
Acrylonitrile butadiene styrene	4.26 - 5.60
Polyester	1.85 - 2.43
Chloroform	<1

NOTE: Accessories containing stainless steel may contain up to 0.4% weight of the CMR substance Cobalt.

For symbols definitions, visit:

<https://www.medcomp.net/patientresources/symbol-glossary>

Medcomp® and Symetrex® are trademarks of Medical Components, Inc. registered in the United States.

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The ID card template includes a header with the URL <http://www.medcomp.net/patientinformation>. It features a graphic of a person with a question mark and a magnifying glass over a box labeled 'ID'. Below this, there are four horizontal lines for text entry, corresponding to the instructions provided. At the bottom, the Medcomp, Inc. logo and address are listed: Medical Components, Inc., 1199 24th St., Hudsonville, MI 49421, U.S.A., P: 215-256-4201, www.medcomp.net.

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The diagram shows a patient ID card with the following fields and icons:

- Top left: A small icon of a person with a question mark.
- Top right: A URL: <http://www.medcomp.net/patientinformation>
- Left side: A vertical column with three icons: a person, a question mark, and a person with a question mark.
- Right side: Four horizontal lines for text entry, corresponding to the fields in the list above.
- Bottom left: A logo for Medical Components, Inc. with address: 1091 Delo Dr., Indianapolis, IN 46240, U.S.A. Phone: 215-256-4201. Website: www.medcomp.net

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Diagram of a patient ID card. The card has a header with a URL: <http://www.medcomp.net/patientinformation>. Below the header are three main sections: 1. A section with a question mark icon and a line for the patient's name. 2. A section with a calendar icon and a line for the implantation date. 3. A section with a person icon and a line for the provider's name and address. At the bottom is a logo for Medial Components, Inc. with the address: 1109 Delp Dr., Hialeahville, PA 19148, U.S.A., Tel: 215-256-4201, and the website: www.medcomp.net.

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Medcomp® ve Symetrex®, Medical Components, Inc. firmasının Amerika Birleşik Devletleri'nde kayıtlı ticari markalarıdır.
Harleysville, PA 19438, U.S.A.
MPS Medical-Product-Service UK Ltd
MPS Medical Product Service GmbH