

T-3[®] CT

SHORT-TERM HEMODIALYSIS CATHETER



FEATURES & BENEFITS

INTERNAL LUMEN DESIGN

Power Injectable Triple



MATERIAL

Polyurethane

TIP DESIGN

Step Tip




FRENCH SIZE

15.5F



CONFIGURATION

Straight



T-3[®] CT

SHORT-TERM HEMODIALYSIS CATHETER

BASIC SET		
MC013501	15.5F X 15CM T-3 [®] CT, STRAIGHT	5/BOX
MC013502	15.5F X 20CM T-3 [®] CT, STRAIGHT	5/BOX
MC013503	15.5F X 24CM T-3 [®] CT, STRAIGHT	5/BOX
MC013504	15.5F X 28CM T-3 [®] CT, STRAIGHT	5/BOX
MC013505	15.5F X 32CM T-3 [®] CT, STRAIGHT	5/BOX

Basic Sets Contents: (1) Catheter (1) Guidewire (1) Guidewire Advancer
(1) Introducer Needle (1) Scalpel (1) Suture (1) Valved Peelable
Introducer (2) Dilator (3) End Cap

IMPORTANT RISK INFORMATION

Indications For Use: The T-3[®] CT Catheter is indicated for short-term use where vascular access is required for less than 14 days for the purpose of hemodialysis. The third internal lumen is indicated for intravenous administration of fluids or medications and power injection of contrast media.

Contraindications: Known or suspected allergies to any of the components of the catheter or the kit. This device is contraindicated for patients exhibiting severe, uncontrolled coagulopathy or thrombocytopenia.

Refer to Instructions for Use for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation



Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201
www.medcomp.net

ISO13485
Certified Quality System

PN2731BSI Rev. B 11/25