




SYMETREX[®]

LONG TERM HEMODIALYSIS CATHETER

Features & Benefits

-  **Double “D”** Internal Lumen Design
-  **Polyurethane** Material
-  **Symmetrical** Tip Design
-  **Straight** Configuration
-  **15.5F** Available

SYMETREX[®]

LONG TERM HEMODIALYSIS CATHETER

SYMETREX[®] WITH SIDHOLES

CATHETER ONLY SET		
MUPS019SE.	15.5F X 19CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MUPS023SE.	15.5F X 23CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MUPS028SE.	15.5F X 28CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MUPS033SE.	15.5F X 33CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
CATHETER SET WITH STYLETS		
MBPS019SE.	15.5F X 19CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MBPS023SE.	15.5F X 23CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MBPS028SE.	15.5F X 28CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MBPS033SE.	15.5F X 33CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
CATHETER SET WITHOUT STYLETS		
MBNS037SE.	15.5F X 37CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MBNS042SE.	15.5F X 42CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX

SYMETREX[®] WITHOUT SIDHOLES

CATHETER ONLY SET		
MUPS019E.	15.5F X 19CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MUPS023E.	15.5F X 23CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MUPS028E.	15.5F X 28CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MUPS033E.	15.5F X 33CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
CATHETER SET WITH STYLETS		
MBPS019E.	15.5F X 19CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MBPS023E.	15.5F X 23CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MBPS028E.	15.5F X 28CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MBPS033E.	15.5F X 33CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
CATHETER SET WITHOUT STYLETS		
MBNS037E.	15.5F X 37CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX
MBNS042E.	15.5F X 42CM SYMETREX [®] CATHETER, STRAIGHT	5/BOX

Catheter Only Set Contents: (1) Catheter (2) Stylet (2) End Cap

Catheter Set Contents: (1) Catheter (2) Stylet* (1) Introducer Needle (1) Guidewire (1) Guidewire Advancer (1) Tunneler (2) Dilator (1) Valved Peelable Introducer (2) End Cap

*Not included in 37CM & 42CM sets

IMPORTANT RISK INFORMATION

Indications for Use: The Symetrex[®] Catheter is indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis.

Contraindications: Known or suspected allergies to any of the components of the catheter or the kit. This device is contraindicated for patients exhibiting severe, uncontrolled coagulopathy or thrombocytopenia. This device is contraindicated for subclavian puncture when ventilator is in use. The presence of other device-related infection, or septicemia is known or suspected. Severe chronic obstructive lung disease is present. Tissue factors in the localized area of device placement will prevent proper device stabilization and/or access. Post irradiation of prospective insertion site.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation

